

The Quickest way to travel



GO GO
ELITE TRAVELLER®

The compact design of the Go-Go Elite Traveller allows you to easily manoeuvre in tight spaces while providing stable outdoor performance. With a 136 kg weight capacity, a charger port in the tiller, one-hand feather-touch disassembly and maximum speeds up to 7 km/h, it's easy to go where you want to go!





Feather-Touch Disassembly

PRIDE®
MOBILITY
LIVE YOUR BEST®

Specifications

3W: SC40E 4W: SC44E

WEIGHT CAPACITY	136 kg
MAXIMUM SPEED ¹	Up to 7 km/h
GROUND CLEARANCE ²	35 mm
TURNING RADIUS	3W: 838 mm 4W: 1120 mm
OVERALL LENGTH ²	3W: 940 mm 4W: 1003 mm
OVERALL WIDTH ²	495 mm
FRONT TIRES	3W: 2" x 8" solid 4W: 2" x 7" solid
REAR TIRES	2.5" x 8" solid
RANGE PER CHARGE ^{1,4}	3W: Up to 11,6 km (12 AH) 4W: Up to 9,7 km (12 AH) 3W: Up to 16,7 km (18 AH) 4W: Up to 15,6 km (18 AH)
TOTAL WEIGHT WITHOUT BATTERIES	3W: 34,2 kg 4W: 36 kg
HEAVIEST PIECE WHEN DISASSEMBLED	3W: 13,4 kg (front section) 4W: 15,4 kg (front section)
STANDARD SEATING	Type: Foldable, compact, vinyl-covered, molded plastic Weight: 10,2 kg Material: Black Vinyl Dimensions: Width: 431 mm Depth: 431 mm
MOTOR OUTPUT	24-volt DC motor, sealed mini-transaxle
BRAKES	Regenerative and electromechanical
BATTERY REQUIREMENTS ^{3,5}	Size: (2) 12 AH Weight: 4 kg each, 9 kg per battery pack Size: (2) 18 AH (optional) Weight: 6,3 kg each, 13,4 kg per battery pack
BATTERY CHARGER	2 amp, external
STANDARD COLOURS	2 sets with every scooter!   Red Blue

Features

- Charger port in tiller
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Includes two sets of easily changeable, red and blue colored panels
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

Accessories

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 1270 mm lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- XLR USB charger

Options

- Additional 12AH battery box
- Additional 18AH battery box

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards.