

# Excellent Outdoor Performance With Luxurious Features



**Luna Victory**  
PRIDE® MOBILITY PRODUCTS





The Pride® Luna Victory® E 3-wheel Scooter gives a powerful performance for both in- and outdoor usage. Luxurious styling and comfort features like a foam delta tiller and non-marking tires make the Pride Luna Victory® E an all in-one-package.



**PRIDE®**  
MOBILITY  
LIVE YOUR BEST®

# Specifications

SC613

|   |   |
|---|---|
| WEIGHT CAPACITY                             | 159 kg  |
| MAXIMUM SPEED <sup>1</sup>                  | Up to 15 km/h   |
| GROUND CLEARANCE <sup>2</sup>               | 95 mm   |
| TURNING RADIUS <sup>2</sup>                 | 1257 mm   |
| OVERALL LENGTH <sup>2</sup>                 | 1397 mm   |
| OVERALL WIDTH <sup>2</sup>                  | 699 mm  |
| FRONT TIRES                                 | 4" x 13" pneumatic  |
| REAR TIRES                                  | 4" x 13" pneumatic  |
| SUSPENSION                                  | Front and rear  |
| DRIVE SYSTEM                                | Rear-wheel drive, sealed transaxle, 24-volt DC Motor  |
| BRAKES                                      | Regenerative and electromechanic  |
| RANGE PER CHARGE <sup>1,4</sup>             | Batteries 60 AH: up to 30 km<br>Batteries 75 AH: up to 40 km  |
| TOTAL WEIGHT WITHOUT BATTERIES <sup>6</sup> | 87,5 kg   |
| HEAVIEST PIECE WHEN DISASSEMBLED            | 35,6 kg (front)   |
| STANDARD SEATING                            | Rotatable and in height adjustable seat provided with seat slider, in height adjustable armrests and adjustable back.<br>Weight: 21 kg<br>Material: Black Vinyl<br>Dimensions: Width: 457 mm<br>Depth: 457 mm |
| BATTERY REQUIREMENTS <sup>4</sup>           | Two 12-volt, deep-cycle<br>Type: 60 AH (standard)<br>73 AH (optional)<br>Weight <sup>5</sup> : 60 AH: 17,5 kg each;<br>75 AH: 24,49 kg each   |
| BATTERY CHARGER                             | 5 amp, off-board  |
| COLOURS                                     |  Red  Yellow                            |

# Features

- Large 13" pneumatic tires for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, LED head and tail lights
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard
- Easy throttle control by switchable push and pull lever
- Steady speed adjustment
- Tiller display operable with just one hand

# Accessories

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 1270 mm lap belt
- Oxygen tank holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Hand brake
- Seat with adjustable lumbar

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Battery weight may vary based on manufacturer.

<sup>5</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.