

# Features, Performance, and Travel Ready



**GO GO**  
ELITE TRAVELLER® PLUS





Feather-touch disassembly

The Go-Go Elite Traveller® Plus offers style and performance for your on-the-go lifestyle. Standard features include a delta tiller with wraparound handles, feather-touch disassembly, black, non-scuffing tires, front basket and an 18 AH battery pack. A 147 kg weight capacity and a maximum speed up to 6.4 km/h complement the versatility of the Go-Go Elite Traveller Plus.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

## Specifications

3W: SC53 4W: SC54

WEIGHT CAPACITY	147 kg
MAXIMUM SPEED <sup>1</sup>	Up to 6,4 km/h
GROUND CLEARANCE <sup>2</sup>	38 mm
TURNING RADIUS <sup>2</sup>	3W: 876 mm   4W: 1162 mm
OVERALL LENGTH <sup>2</sup>	3W: 1016 mm   4W: 1023 mm
OVERALL WIDTH <sup>2</sup>	540 mm
FRONT TIRES	2" x 8" solid
REAR TIRES	3" x 9" solid
RANGE PER CHARGE <sup>1,4</sup>	Up to 20,25 km
TOTAL WEIGHT WITHOUT BATTERIES	3W: 39,2 kg   4W: 39,5 kg
HEAVIEST PIECE WHEN DISASSEMBLED	3W: 14,7 kg (front section) 4W: 15,9 kg (front section)
STANDARD SEATING	Type: Foldable, compact, vinyl-covered, molded plastic Weight: 10,2 kg Material: Black Vinyl Dimensions: Width: 431 mm Depth: 431 mm
MOTOR OUTPUT	24-volt DC motor, sealed mini-transaxle
BRAKES	Regenerative and electromechanical
BATTERY REQUIREMENTS <sup>3,5</sup>	Two 12-volt, deep cycle (standard) Size: 2 x 18 AH (standard) Weight: 18 AH: 13.2 kg (battery pack)
BATTERY CHARGER	2 amp, off-board
STANDARD COLOURS	2 sets with every scooter!   Red Blue

## Features

- Feather-touch disassembly permits simple frame separation with only one hand and easily
- Delta tiller with wraparound handles
- Pride's exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- Includes two sets of easily changeable, red and blue colored panels
- Front frame-mounted seat post offers maximum stability
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Disassembles into five lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Front tiller-mounted basket
- Removable footboard basket for additional storage (4-wheel version only)
- Convenient off-board dual voltage charger can charge battery pack on- or off-board

## Accessories

- Additional 18AH battery box
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 1270 mm lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
- <sup>3</sup> AGM or gel cell type required.
- <sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>5</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards.