

# The ultimate High-performance scooter





**GO-GO**<sup>®</sup>  
ELITE TRAVELLER<sup>®</sup> *SPORT*

The stylish, high-performing Go-Go Elite Traveller<sup>®</sup> Sport features an innovative design with a full complement of unique standard features. A charger port conveniently located in the tiller, front and rear LED lighting, 9" wheels at front and rear, an additional 50 mm deck space and an increased weight capacity of 147 kg are just a few features that make the Go-Go Elite Traveller Sport an exceptional scooter choice.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

## Specifications

SC74

WEIGHT CAPACITY	147 kg
MAXIMUM SPEED <sup>1</sup>	Up to 6,4 km/h
GROUND CLEARANCE <sup>2</sup>	38 mm
TURNING RADIUS <sup>2</sup>	1283 mm
OVERALL LENGTH <sup>2</sup>	1092 mm
OVERALL WIDTH <sup>2</sup>	540 mm
FRONT TIRES	9" solid
REAR TIRES	9" solid
RANGE PER CHARGE <sup>1,4</sup>	Up to 20,25 km
TOTAL WEIGHT WITHOUT BATTERIES	44 kg
HEAVIEST PIECE WHEN DISASSEMBLED	20 kg (front)
STANDARD SEATING	Type: Compact foldable molded plastic Weight: 10 kg Material: Black vinyl Dimensions: 430 mm overall width 430 mm effective depth
MOTOR OUTPUT	Rear-wheel drive, sealed transaxle, 24-volt DC motor
BRAKES	Regenerative and electromechanical
BATTERY REQUIREMENTS <sup>3,5</sup>	Two 12-volt, deep cycle (standard) Size: 2 x 18 AH (standard) Weight: 18 AH: 13 kg battery pack
BATTERY CHARGER	2-amp, off-board
STANDARD COLOURS	2 sets of colors with every scooter!   Red Blue

## Features

- Front and rear LED lighting
- Feather-touch disassembly with one hand
- Additional 2 inches of deck space for added leg room
- Exclusive black, low-profile, non-scuffing tires with new tri-spoke mags
- Delta tiller with wraparound handles
- Auto-connecting front to rear cable
- Includes 2 sets of easily changeable colored panels - Red & Blue
- Front frame-mounted seat post offers maximum stability
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Front tiller-mounted basket
- Charger port conveniently located in the tiller
- Suspension seat post

## Accessories

- Additional 18AH battery box
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 1270 mm lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Includes standard seat.

<sup>5</sup> Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards.