

Simply Fold and Travel




Folded Size
787 x 406 x 600 mm

Experience the beauty of lightweight travel with the iGo+. This compact power chair folds in just a few simple steps for effortless transport. The iGo+ makes traveling a breeze with its convenient storage.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

WEIGHT CAPACITY	113 kg
MAXIMUM SPEED ¹	Up to 6 km/h
GROUND CLEARANCE ³	79 mm
TURNING RADIUS ³	616 mm
OVERALL LENGTH ³	1016 mm
OVERALL WIDTH ³	597 mm
SEAT-TO-GROUND HEIGHT	483 mm
SEAT-TO-FOOT PLATFORM HEIGHT	368 mm
CASTERS	8" solid
DRIVE TIRES	12" pneumatic drive
SUSPENSION	Limited
RANGE PER CHARGE ^{1,2}	Up to 15 km
TOTAL WEIGHT WITHOUT BATTERIES ^{3,5}	24 kg
HEAVIEST PIECE WHEN DISASSEMBLED ³	N/A
STANDARD SEATING	Material: Gray, Foam Dimensions: 445 x 419 mm
DRIVETRAIN	2x 180W brushed motors
DUAL BRAKING SYSTEM	Regenerative and electromechanical
BATTERY REQUIREMENTS ³	24V 18Ah Lithium ion
BATTERY WEIGHT ⁴	3,3 kg
BATTERY CHARGER	2-amp, off-board
COLOR	 Silver

Features

- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Includes convenient, under seat mesh storage bag, mesh cup holder and 1016 mm lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort

Accessories

- XLR USB charger



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

⁵ Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Pride Mobility Products Europe BV
De Zwaan 3
1601 MS Enkhuijzen
The Netherlands

Tel: +31 (0)85-0408630 • Fax: +31 (0)85-0430039
quantumrehab.eu • info@prideeurope.nl