

Elevate your life to the next level





Jazzy air²

At the touch of a switch, the Jazzy Air[®] 2 elevates 300 mm in just 11 seconds. With a sleek, more modern look, this power chair safely drives at speeds up to 6,5 km/h while elevated, allowing users to enjoy face-to-face social engagement. With Mid-Wheel 6[®] Drive and Active-Trac[®] Suspension, the Jazzy Air[®] 2 offers enhanced performance and stability both indoors and out.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

WEIGHT CAPACITY	136 kg
MAXIMUM SPEED ¹	Up to 6,4 km/h elevated and non-elevated position
GROUND CLEARANCE ²	25 mm
SEAT TO FLOOR	539 - 838 mm
TURNING RADIUS ²	527 mm
OVERALL SIZE ^{2,6}	Length: 1065 mm Width: 589 mm
DRIVE CONTROLLER	PG VR2
DRIVE WHEELS	10" flat free
CASTER WHEELS	6" flat free
SUSPENSION	Active-Trac® with Mid-Wheel 6* drive Technology
RANGE ^{1,4}	Up to 27 km (U1) Up to 30,9 km (40 Ah)
COMPONENT WEIGHTS ⁵	Base: 104,3 kg ⁷ with U1 Accus 111,1 kg ⁷ with 40 Ah Accus U1 Accus: 11,1 kg 40 Ah Accus: 14,7 kg
DRIVETRAIN	Two motor, Mid-Wheel 6* drive
BRAKES	Intelligent Braking (electronic regenerative, disc park brakes)
BATTERIES ³	(2) 12 V, U1 (2) 40 Ah
BATTERY CHARGER	5 amp, external
WARRANTY	5 year limited (frame) 13 month limited (batteries) 13 month limited (electronics) 13 month (drive motors)
COLORS	 <p>Citrine Yellow Sapphire Blue Pink Topaz Onyx Black Ruby Red Pearl White Silver</p>
MATTE COLORS	 <p>Black Pearl Tanzanite Garnet Red</p>

Features

- High-back or solid seat pan, reclining sport seat. Seat features 50 mm of depth adjustment. Choice of 400x400-450 mm, 450x450-500mm or 500x500-560mm
- Angle, height and depth adjustable armrests
Height adjustable headrest
- Newly re-designed lift features:
 - 50 mm lower seat-to-floor height enables user to sit comfortably at standard table or desk
 - Raises and lowers in 11 seconds by using the Air switch or VR2 joystick
- Drives safely at 6,4 km/h in the elevated position comfort in both the elevated and non-elevated positions
- Bright LED marker lights offers visibility in all environments
- Active-Trac® suspension with Mid-Wheel 6* Technology allows for maximum stability in the elevated and non-elevated positions.
- Single-sided caster forks for enhanced maneuverability
- Under seat front storage compartment keeps your valuables accessible

Accessories

- Alternative size lap belt: 1270 mm, 1780 mm, 2300 mm or 2280 mm
- Cane holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB Charger

Options

- 40 amp AGM batteries
- Swing-away joystick
- Stealth® cushions

1) Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%. Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.

2) AGM or Gel-Cell type required.

3) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

4) Battery weight may vary based on manufacturer.

5) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

7) Weight includes base, seat and batteries.

NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, and EN12184 standards. All specifications subject to change without notice.