

The Ultimate Compact and Foldable travel scooter



QUEST






Ready for transport and easily folded; the 4-wheel Quest Folding Scooter is perfect for the active lifestyle. The Quest provides a comfortable and convenient ride offering great performance in- and outdoors at a maximum speed of 6,4 km/h and a range per charge of up to 20 km. It's a super portable space saver that adjusts to fit perfectly inside any closet, corner or vehicle trunk. The Quest folding scooter is sleek, stylish, and extremely easy to operate, making it the perfect travel mobility scooter choice.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

QUEST

Specifications

S19

WEIGHT CAPACITY	113,4 kg
MAXIMUM SPEED ¹	Up to 6 km/h
GROUND CLEARANCE ²	58 mm
TURNING RADIUS ²	813 mm
OVERALL LENGTH ²	940 mm
OVERALL WIDTH ²	483 mm
SEAT-TO-GROUND HEIGHT	520 mm
SEAT-TO-DECK HEIGHT	330 mm
FRONT TIRES ²	1.5" x 7" solid
REAR TIRES ²	2" x 8" solid
SUSPENSION	N/A
RANGE PER CHARGE ^{1,3}	Up to 20 km
TOTAL WEIGHT WITHOUT BATTERIES	23,5 kg
HEAVIEST PIECE WHEN DISASSEMBLED	N/A (unit doesn't disassemble)
STANDARD SEATING	Type: Foldable, compact, molded plastic Material: Gray Vinyl Dimensions: 406 x 330 mm
DRIVE SYSTEM	4-Pole 270w
DUAL BRAKING SYSTEM	Regenerative and electromechanical
BATTERY REQUIREMENTS ⁴	Size: (2) 12 AH (sealed lead acid) standard Weight: 4 kg each (2 required), or Size: (1) 11.5 AH (lithium-ion optional) Weight: 2,3 kg each (1 required)
BATTERY CHARGER	2-amp, off-board (sealed lead acid) or 2.5-amp, off-board (lithium-ion)
COLORS	 White  Black  Red

Features

- Unique folding design
- Folds in just three easy steps
- Height adjustable tiller
- Up to 20 km range per charge
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down
- Articulated front end adds greater stability
- Front and rear LED lights
- Optional armrests
- Folded: 762mm, 483mm, 457mm (HxWxD)

To fold:



Fold Seat



Pull Up Lever to Fold



Fold and Lock Tiller

To unfold:



Pull Lever



Unfold Scooter



Unfold Seat; Lock in Place



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Quest folding scooter and its lithium-ion battery are permitted.

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- ³ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- ⁴ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.