

Excellent Outdoor Performance with Luxurious Features



Luna Victory
PRIDE® MOBILITY PRODUCTS





The Pride Luna Victory E 3-wheel scooter gives a powerful performance for both in- and outdoor usage. Luxurious styling and comfort features like a foam delta tiller and non-marking tires make the Pride Luna Victory E an all in-one-package.



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Specifications

SC613

WEIGHT CAPACITY	159 kg
MAXIMUM SPEED ¹	Up to 12,6 km/h
GROUND CLEARANCE ²	95 mm
TURNING RADIUS ²	1257 mm
OVERALL LENGTH ²	1403 mm
OVERALL WIDTH ²	686 mm
FRONT TIRES	13" pneumatic
REAR TIRES	13" pneumatic
SUSPENSION	Front and rear
DRIVE SYSTEM	Rear-wheel drive, sealed transaxle, 24-volt DC Motor
BRAKES	Regenerative and electromechanic
RANGE PER CHARGE ^{1,4}	Batteries 60 AH: up to 30 km Batteries 75 AH: up to 40 km
TOTAL WEIGHT WITHOUT BATTERIES ⁶	88 kg
STANDARD SEATING	Type: High Back Lumbar Weight: 23 kg Material: Black Vinyl Width: 457 mm Depth: 483 mm
BATTERY REQUIREMENTS ⁴	Two 12-volt, deep-cycle Size: 70 Ah or 55 Ah Weight ⁵ : 60 AH: 17,5 kg each; 75 AH: 24,49 kg each
BATTERY CHARGER	5 amp, off-board
COLOURS	  Red Yellow

Features

- Large 13" pneumatic tires for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, LED head and tail lights
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard
- Easy throttle control by switchable push and pull lever
- Steady speed adjustment
- Tiller display operable with just one hand

Accessories

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 1270 mm lap belt
- Oxygen tank holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Hand brake
- Seat with adjustable lumbar

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.

³ AGM or gel cell type required.

⁴ Battery weight may vary based on manufacturer.

⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.