

Outstanding Outdoor Performances with luxurious features




VICTORY *XL 130*

The **Victory XL 130** features a powerful drivetrain and large 13" wheels to handle rugged outdoor terrain. Luxury touches like a wraparound delta tiller, low-profile nonmarking pneumatic tires and a LED headlight blend seamlessly give the Victory XL 130 it's sleek, sporty style.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

Model number	SC713
Weight capacity ⁵	181 kg
Maximum speed ¹	Up to 16 km/h
Ground clearance ²	90 mm
Turning radius ²	1695 mm
Length ²	1353 mm
Width ²	641 mm
Tires	13" pneumatic (Front and rear)
Suspension	Front and rear
Range per charge ^{1,3}	38 km
Total weight without batteries ⁵	78,5 kg
Standard seating	Type: High Back Material: Black Vinyl Width: 450 mm Depth: 450 mm Weight: 20,75 kg
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system	Regenerative and electromechanical
Battery requirements ³	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 55 or 75 Ah
Battery weight ⁴	55 Ah: 17 kg 75 Ah: 24,4 kg
Battery charger	5 or 8 Amp, off-board
Colors	 <p>Candy Apple Red Onyx Black Mercury Silver Viper Blue</p>

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3. AGM or gel cell type required.
4. Battery weight may vary based on manufacturer.
5. Including standard seat

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Features

- Large 13" pneumatic tires for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable LED headlight
- Delta tiller with wraparound handles
- Front basket standard
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard

Accessories

- Black weather cover
- Rearview mirror
- Cup holder
- Lap belt
- Safety flag
- Saddle bag
- Cell phone holder
- XLR USB charger
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane holder, single or double

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressure, ambient temperature, and the gradient and surface of the road or pavement.