

PRODUCT INSERT

POWER LIFT RECLINER WEIGHT CAPACITIES

Your power lift recliner was developed for your convenience. See the below specification table for information on weight capacity. Refer to the owner's manual for any additional information concerning your power lift recliner.

| MEDICAL POWER LIFT RECLINER WEIGHT CAPACITY | |
|---|------------------------------------|
| MODEL | WEIGHT CAPACITY |
| C1 and C11 | 113 kg/17 stone 12 lbs. (250 lbs.) |
| 101 and 107 | 147 kg/23 stone 3 lbs. (325 lbs.) |
| D30, 660 and 670 | 170 kg/26 stone 11 lbs. (375 lbs.) |



MANDATORY! Stay within the specified weight capacity of your power lift recliner. Exceeding the weight capacity voids your power lift recliner warranty.

NOTE: *Power Lift Recliner model numbers may contain a prefix or suffix (i.e., prefixes LC, NM or SR; suffixes S, M, L, PW, etc.). When using the table above, look for only the numeric portion of your specific model.*

