

Economy class price First class features









GO GO[®]
ULTRA X

The Go-Go[®] Ultra X is loaded with features such as one-hand disassembly and a convenient drop-in battery box for worry-free travel. With a 117 kg weight capacity, a maximum speed up to 6,4 km/h, a per-charge range up to 17,5 km on the 3-wheel and 18,4 km on the 4-wheel, Go-Go[®] Ultra X is an exceptional value.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Specifications

3W: 40X 4W: 44X

Weight capacity	118 kg
Max. speed ¹	Up to 6,4 km/h
Ground clearance ²	32 mm
Turning radius ²	3W: 781 mm 4W: 1385 mm
Length ²	3W: 940 mm 4W: 1120 mm
Width ²	495 mm
Front tires	3W: 8" solid 4W: 7" solid
Rear tires	8" solid
Max. range per charge ^{1,4}	3W: Up to 11 km 4W: Up to 11,5 km
Total weight without batteries ⁵	3W: 34,5 kg 4W: 36,3 kg
Standaard seat ²	Type: Folded Material: Black Vinyl Width: 432 mm Depth: 432 mm Weight: 10,4 kg
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Batteries	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 12, 18 or 22 Ah
Standard colors	  Red Blue
Optional colors	    Orion Yellow Onyx Black Pearl Pink Sunburst Orange

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Features

- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Pride's exclusive black, non-scuffing tires

Accessories

- Cane/crutch holder (single, double)
- Cup holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger