

# Supreme performance, exhilarating style and **Excellent Value**





The sleek, sporty Victory® 10 DX Scooter delivers smooth performance and all-new features. With a full complement of advanced standard features including rugged suspension, a high-visibility automatic front and rear lighting system, feather-touch disassembly, an LED curb light, and a wraparound delta tiller, the Victory® 10 DX is a greater value than ever before.



## Specifications

3W: S610DX 4W: SC710DX

Weight capacity <sup>7</sup>	159 kg
Maximum speed <sup>1</sup>	Up to 10 km/h
Ground clearance <sup>2</sup>	63 mm
Turning radius <sup>2</sup>	1067 mm
Length <sup>2</sup>	3W: 1187 mm   4W: 1194 mm
Width <sup>2</sup>	562 mm
Tires	10" solid (front and rear)
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DX motor
Dual braking system	Regenerative and electromechanical
Range per charge <sup>1,4</sup>	Up to 25 km
Total weight without batteries <sup>5</sup>	3W: 66,9 kg   4W: 84,8 kg
Standard seating	Type: CRS Junior Weight: 19 kg Material: Black Vinyl Width: 445 mm Depth: 432 mm
Batteries <sup>3,6</sup>	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: U1 (31, 32 or 36 Ah) Weight: 10 kg
Battery charger	5 Amp, off-board
Colors	  Red Blue

1. Varies with user weight, terrain type, battery amp-hour rating (Ah), battery charge, battery condition, tire type and tire condition. This specification can be subject to a variance of +10%, -5%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. AGM or Gel-Cell type required.
4. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Includes standard seat.
6. Battery weight may vary based on manufacturer.
7. All measurements taken using standard seat. Measurements may vary with use of different seating styles.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Features

- Limited front & full rear suspension
- High visibility automatic front & rear lighting system
- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Long-lasting, high intensity LED headlight, tail light & indicators
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- New backlit battery meter for increased accuracy & visibility in any setting
- CRS high-back seat with headrest
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port for standard dual voltage off-board charger
- Rearview mirror standard
- Seat post suspension standard

## Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover