

4-wheel stability with 3-wheel maneuverability



ZERO TURN 8

(also known as the Jazzy Zero Turn)



Get a 4-wheel scooter with a 3-wheel turning radius with Zero Turn 8 and Pride's patented, intelligent turning, iTurn Technology. Navigate tight corners and small spaces effortlessly with a 970 mm turning radius. Illuminate your ride with bright LED lighting. Dualmotors and CTS suspension ensure a smooth and comfortable drive.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

Jazzy ZTS

Weight capacity	147 kg
Max. speed ¹	Up to 6,4 km/h
Ground clearance ²	33 mm
Turning radius ²	948 mm
Length ²	1162 mm
Width ²	525 mm
Tires	Front: 7" solid Rear: 9" solid
Range per charge ^{1,4}	Up to 17,7 km at 147 kg Up to 21,1 km at 90,7 kg
Total weight without batteries ⁵	54 kg
Standard seat ²	Type: Foldable, compact, molded plastic Weight: 12 kg Material: Black vinyl Width: 432 mm Depth: 432 mm
Drive system	Rear-wheel drive, Two motor, 24-volt DC motor
Battery requirements	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 18 or 22 Ah Weight: 5,25 kg each
Colors	

1. Varies with user weight, terrain type, battery amp-hour rating (Ah), battery charge, battery condition, tire type and tire condition. This specification can be subject to a variance of +10%, -5%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. AGM or Gel-Cell type required. See III. "Batteries and Charging".
4. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Includes standard seat.
6. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Features

- The first 4-wheel indoor scooter
- Charging port built into the tiller for easy charging of the scooter
- Contemporary rim design makes the Zero Turn 8 look as good as it drives
- Dual motors for added traction in rough terrain
- Convenient under seat storage keeps your personal items close at hand
- USB charger built into the tiller for convenient charging of smart phones and portable devices

Accessories

- Rear basket
- Oxygen holder
- Cell phone holder
- Cup holder
- Crutch holder
- Lap belt
- Weather cover
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder