

Features, Performance and Travel Ready





GO GO
ELITE TRAVELLER® PLUS

The Go-Go Elite Traveller Plus offers style and performance for your on-the-go lifestyle. Standard features include a delta tiller with wraparound handles, feather-touch disassembly, black, non-scuffing tires and a standard front basket. The 147 kg weight capacity and a maximum speed up to 6,4 km/h complement the versatility of the Go-Go Elite Traveller Plus.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Specifications

3W: SC53 4W: SC54

Weight capacity	147 kg
Max. speed ¹	Up to 6,4 km/h
Ground clearance ²	38 mm
Turning radius ²	3W: 876 mm 4W: 1162 mm
Length ²	3W: 1016 mm 4W: 1023 mm
Width ²	540 mm
Tires	3W: 8" solid 4W: 9" solid
Max. range per charge ^{1,4}	Up to 20,25 km
Total weight without batteries ⁵	3W: 39,25 kg 4W: 44 kg
Standaard seat ²	Type: Folded Material: Black Vinyl Width: 432 mm Depth: 432 mm Weight: 10,2 kg
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Batteries ⁵	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 12, 18 or 22 Ah Weight: 5,25 kg each
Standard colors	  Red Blue

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication, we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Features

- Feather-touch disassembly permits simple frame separation with only one hand and easily
- Delta tiller with wraparound handles
- Pride's exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- Includes two sets of easily changeable, red and blue colored panels
- Front frame-mounted seat post offers maximum stability
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Disassembles into five lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Front tiller-mounted basket
- Removable footboard basket for additional storage (4-wheel version only)
- Convenient off-board dual voltage charger can charge battery pack on- or off-board

Accessories

- Cane/crutch holder (single, double)
- Cup holder
- Oxygen tank holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger