

The ultimate High-Performance scooter



GO GO
ELITE TRAVELLER[®] *SPORT*



The stylish, high-performing Go-Go Elite Traveller[®] Sport features an innovative design with a full complement of unique standard features. A charger port conveniently located in the tiller, front and rear LED lighting, 9" wheels at front and rear, an additional 50 mm deck space and an increased weight capacity of 147 kg are just a few features that make the Go-Go Elite Traveller[®] Sport an exceptional scooter choice.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Specifications

SC74

Features

Weight capacity	147 kg
Maximum speed ¹	Up to 6,4 km/h
Ground clearance ²	38 mm
Turning radius ²	1283 mm
Length ²	1092 mm
Width ²	540 mm
Tires	9" solid
Range per charge ^{1,4}	Up to 19 km
Total weight without batteries	44 kg
Heaviest piece when disassembled	Front section: 20 kg
Standard seating	Type: Folded Material: Black Vinyl Width: 432 mm Depth: 432 mm Weight: 10,2 kg
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery requirements ^{3,5}	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 12, 18 or 22 Ah Weight: 5,25 kg each
Colors	  Red Blue

- Front and rear LED lighting
- Feather-touch disassembly with one hand
- Additional 2 inches of deck space for added leg room
- Exclusive black, low-profile, non-scuffing tires with new tri-spoke mags
- Delta tiller with wraparound handles
- Auto-connecting front to rear cable
- Includes 2 sets of easily changeable colored panels – Red & Blue
- Front frame-mounted seat post offers maximum stability
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Front tiller-mounted basket
- Charger port conveniently located in the tiller
- Suspension seat post

Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
3. AGM or gel cell type required.
4. Includes standard seat.
5. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.