

Durable, Lightweight and Portable




Folded Size
630mm x 355mm x 770mm

The i-Go from Pride Mobility features an advanced folding technology, enabling it to be transported with ease. The i-Go features front suspension, a durable seating system, under seat storage and much more. Designed to fit inside any small space, the i-Go is the perfect choice for the active individual.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

Weight Capacity	100 kg
Maximum Speed ²	Up to 4,5 km/h
Turning Radius ¹	900 mm
Overall Size ¹	900 mm x 600 mm x 880 mm
Seat Height	500 mm
Seat Width	400 mm
Seat Depth	410 mm
Tires	Front: 8,5" solid Rear: 5,7" solid
Standard Controls	VR2 controller
Range per Charge	Up to 10 km
Total Weight	19,8 kg without batteries
Battery Requirements ¹	Quantity: 1 + 1 spare Type: Lithium battery
Battery Weight ^{1,3}	1,4 kg
Colour	 Silver

¹ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %.

² Specification varies with user weight, terrain, battery type, battery condition and tire pressure.

³ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Features

- Quickly folds in 5 easy steps
- Unique seating for added comfort
- Front suspension
- Large, secure storage compartment
- Lightweight Lithium-ion batteries
- Easy-to-use VR2 controller

How to fold



1 Switch off and remove joystick



2 Place joystick in backrest pocket and remove cushion



3 Fold backrest down, then pull both arm rests anticlockwise toward the back



4 Pull the two black toggles (located on the back) together



5 Push down on the seat frame until unit fully flattens



6 The i-Go is now ready to go!