

Advanced engineering for **Powerful Performance** and exceptional style



VICTORY XL 140

Bristling with high-performance features like full suspension, low-profile tires and a hydraulic sealed brake system, the Victory XL 140 is a blend of power and precision. The Victory XL 140 features sleek style and a wealth of standard luxury touches like a wraparound delta tiller, new double A-arm suspension, new terrain gripping pneumatic tires and LED lighting, making it the ultimate outdoor scooter.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

SC714

Features

| | |
|---|---|
| Weight capacity ⁵ | 159 kg |
| Maximum speed ¹ | Up to 15 km/h |
| Ground clearance ² | 64 mm |
| Turning radius ² | 1803 mm |
| Length ² | 1422 mm |
| Width ² | 692 mm |
| Tires | 13" pneumatic (Front and rear) |
| Suspension | Front and rear |
| Range per charge ^{1,4} | Up to 48 km |
| Total weight without batteries ⁶ | 123 kg |
| Standard seating | Type: High Back Lumbar Weight: 23 kg Material: Black Vinyl Width: 457 mm Depth: 483 mm |
| Drive system | Rear-wheel drive, sealed transaxle, 24-volt DC motor |
| Dual braking system | Regenerative and electromechanical |
| Battery requirements ³ | Type: Two 12-volt, sealed lead-acid, non-spillable deep-cycle Size: 75 Ah |
| Battery charger | 8 amp, off-board |
| Colors |  Red Blue Black Silver |

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
- AGM or gel cell type required.
- Incl. standard seat

NOTE: All specifications subject to change without notice.

- New Double A-Arm suspension for optimum handling
- Sealed hydraulic brakes
- Large 13" terrain gripping pneumatic tires for excellent outdoor performance
- Full lighting package including directional signals and lower headlight
- Delta tiller with wraparound handles
- Easily accessible tie down points (for transport of unoccupied scooter)
- Rearview mirror standard

Accessories

- Front shield
- Pillow top
- Double crutch holder
- Oxygen tank holder
- Rear basket
- Additional rearview mirror

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressure, ambient temperature, and the gradient and surface of the road or pavement.