



APEX
ALUMA-LITE



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

Experience the convenience of a lightweight, compact travel scooter. Easy to put into the trunk, with its heaviest part weighing only 11 kg.



Maximum safety



Easy-to-use cockpit



Easy to disassemble

Specifications

Weight capacity	120 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 18,5 km
Turning radius	1385 mm
Ground clearance	32 mm
Length	1022 mm
Width	490 mm
Tires	8" solid, front and rear
Drive	Rear-wheel drive, single motor, 24-volt DC motor
Weight without batteries	28,8 kg
Batteries	Type: Single 24-volt Lithium ion Size: 12 Ah
More?	Download our specifications sheet on the website

Features

- Disassembles into 5 pieces
- Airline friendly batteries
- Lightweight aluminum frame
- 8" alloy style wheels
- Easy-to-grip tiller adjustment knob

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors

