

GO GO

ELITE TRAVELLER *SPORT*



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

The stylish Go Go® Elite Traveler Sport offers excellent performance, an innovative design and unique features.



Delta tiller with wraparound handles



9" solid tires



Disassembles into 5 pieces

Specifications

Weight capacity	147 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 19 km
Turning radius	1283 mm
Ground clearance	38 mm
Length	1092 mm
Width	540 mm
Tires	9" solid, front and rear
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Weight without batteries	44 kg
Batteries	Type: Two 12-volt, sealed lead-acid, non-spillable deep-cycle Size: 18 Ah
More?	Download our specifications sheet on the website

Features

- Front and rear LED lighting
- Exclusive black, low-profile, non-suffing tires with new tri-spoke mags
- Delta tiller with wraparound handles
- Modular design for easy serviceability
- Suspension seat post

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors

