

**GO GO** *M<sup>2</sup>*



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039  
pridemobility.eu . info@prideeurope.nl

With its advanced magnesium frame, the Go Go MG redefines what is possible in electric mobility. This unique frame not only provides exceptional stability, but also contributes significantly to minimizing overall weight.



Sleek Design



8" pneumatic tires



Comfort seating

## Specifications

Weight capacity	125 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 10 km
Ground clearance	75 mm
Length	1041 mm
Width	429 mm
Tires	8" pneumatic, front and rear
Drive	Rear-wheel drive
Weight without batteries	14,8 kg
Batteries	Type: Lithium ion Size: 10 Ah
More?	Download our specifications sheet on the website

## Features

- Advanced Magnesium Frame (Only 14.8kg)
- Folded in seconds, perfect for travel and leisure
- 1x 10 Ah lithium battery (suitable for taking on the plane)
- Sleek design
- Available in black and grey

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Available colors

