

GO GO[®]

ULTRA X



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

The Go Go® Ultra X offers the performance you would expect from the number 1 in travel mobility. This makes it the ultimate travel scooter, ready for all your adventures!



7" solid tires on the front



built of strong material



Easy to disassemble

Specifications

Weight capacity	118 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 11,5 km
Turning radius	1120 mm
Ground clearance	32 mm
Length	1003 mm
Width	495 mm
Tires	7" solid front, 8" solid rear
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Weight without batteries	36,3 kg
Batteries	Type: Two 12-volt, deep-cycle Size: 12 Ah
More?	Download our specifications sheet on the website

Features

- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



Red

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Pride Mobility Products Europe
De Zwaan 3
1601 MS Enkhuizen
The Netherlands

Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
pridemobility.eu · info@prideeurope.nl



EU_SSP0011/RevB/01MAR24