

# The Quickest way to travel





The compact design of the Go Go Elite Traveller allows you to easily manoeuvre in tight spaces while providing stable outdoor performance. With a 136 kg weight capacity, a charger port in the tiller, one-hand feather-touch disassembly and maximum speeds up to 6,4 km/h, it's easy to go where you want to go!

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

## Specifications

3W: SC40E 4W: SC44E

Weight capacity	136 kg
Max. speed <sup>1</sup>	Up to 6,4 km/h
Ground clearance <sup>3</sup>	35 mm
Turning radius <sup>3</sup>	3W: 838 mm   4W: 1120 mm
Length <sup>3</sup>	3W: 940 mm   4W: 1003 mm
Width <sup>3</sup>	495 mm
Front tires	3W: 8" solid   4W: 7" solid
Rear tires	8" solid
Max. range per charge <sup>1,2</sup>	3W: Up to 16,7 km 4W: Up to 15,6 km
Total weight without batteries <sup>4</sup>	3W: 34,2 kg 4W: 36 kg
Standard seating	Type: Folded Material: Black Vinyl Width: 432 mm Depth: 432 mm Weight: 10,2 kg
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Batteries <sup>4</sup>	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 12 or 18 Ah
Battery weight (each)	12 Ah: 8,8 kg 18 Ah: 13,5 kg
Colors	  Red Blue

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
4. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication, we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

## Features

- Delta tiller with charger port to conveniently charge your scooter
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces
- Front frame-mounted seat post offers maximum stability
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

## Accessories

- Cane/crutch holder (single, double)
- Cup holder
- Double crutch holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger