

Rugged with feather-touch disassembly



REVO™ 2.0





Experience the **Revo 2.0**, a rugged scooter that combines an array of features into one great package. Tough and built to last, the Revo 2.0 offers rugged dependability you expect from luxury mid-size scooters. Comfort-Trac Suspension (CTS) ensures a smooth and comfortable ride over varied terrain. Experience convenient, feather-touch disassembly for portability that is comparable to our Travel Mobility line. Enjoy standard under seat storage and speeds up to 8,4 km/h. Ride on with the Revo 2.0!

PRIDE®
MOBILITY
LIVE YOUR BEST®



Specifications

Product number	3W: S66 4W: S67
Weight capacity	181 kg
Maximum speed ¹	Up to 8,4 km/h
Ground clearance ²	32 mm
Turning radius ²	3W: 1016 mm 4W: 1323 mm
Length ²	3W: 1167 mm 4W: 1194 mm
Width ²	525 mm
Front tires	3W: 8" solid 4W: 9" pneumatic
Rear tires	3W: 9" solid 4W: 9" pneumatic
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system	Elektromechanic
Range per charge ^{1,4}	Up to 19,5 km at 181 kg Up to 26 km at 91 kg
Total weight without batteries ⁵	3W: 54,25 kg 4W: 59,75 kg
Standard seating	Type: Molded Plastic Material: Black Vinyl Width: 508 mm; 432 mm Depth: 457 mm; 432 mm Weight: 12,25 kg
Batteries ³	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 30 Ah
Battery weight ⁴	14 kg
Battery charger	3.5 amp, off-board
Colors	  Grey Street True Blue

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
4. Battery weight may vary based on manufacturer.
5. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Features

- Comfort-Trac front and rear suspension
- Durable and stylish rear bumper
- Easy accessible tie down points (for transport of unoccupied scooter)
- Easy adjustable tiller
- Easily charge the batteries on or off the unit
- Seat designed for extreme comfort
- Standard front basket
- Standard under seat storage
- Ultra heavy-duty drivetrain with
- User-friendly console
- USB mobile device charging port

Accessories

- Black weather cover
- Rearview mirror
- Cup holder
- Lap belt
- Safety flag
- Saddle bag
- Cell phone holder
- XLR USB charger
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane holder, single or double

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.