



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039  
pridemobility.eu . info@prideeurope.nl

The Colt Sport combines sporty style with advanced features, including a DX Tiller, LED headlights and an advanced suspension seatpost, along with an all-new adjustable memory foam seat.



High-Intensity LED lights



10" solid tires



Easy-to-use cockpit

## Specifications

Weight capacity	159 kg
Speed	Up to 12,9 km/h
Range per charge	Up to 26,8 km
Turning radius	1334 mm
Ground clearance	57 mm
Length	1194 mm
Width	657 mm
Tires	10" solid, front and rear
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Weight without batteries	68 kg
Batteries	Type: Two 12-volt, sealed lead-acid, non-spillable, deep-cycle Size: 40 Ah
More?	Download our specifications sheet on the website

## Features

- Front and rear suspension for optimal outdoor driving on various terrains
- Memory foam adjustable seat
- High-Intensity LED headlight offers optimal pathway illumination including directional signals leg room.
- DX Tiller
- 10" Alloy Style solid tires

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Available colors

