

# GO CHAIR®



**PRIDE®**  
MOBILITY  
LIVE YOUR BEST®

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039  
pridemobility.eu . info@prideeurope.nl

Small, light and portable

When it comes to electric wheelchairs, we changed up the game with the Go Chair.



6 colors available!



Disassembles into 5 pieces



Fits in the car trunk

## Specifications

Weight capacity	136 kg
Speed	Up to 6 km/h
Range per charge	Up to 21,3 km
Turning radius	648 mm
Ground clearance	45 mm
Length	725 - 842 mm
Width	575 mm
Drive wheels	9" solid
Caster wheels	6" solid
Weight without batteries	44,3 kg
Batteries	Type: Two 12 volt, sealed lead-acid, non-spillable deep-cycle Size: 18 Ah
More?	Download our specifications sheet on the website

## Features

- Choice of 6 colours
- Increased 136 kg weight capacity
- Max speed of up to 6 km/h
- Dual swing-away, under-seat storage bins
- One hand feather-touch disassembly

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Available colors

