

# Jazzy<sup>®</sup> EVO 613Li



**PRIDE<sup>®</sup>**  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039  
pridemobility.eu . info@prideeurope.nl

The Jazzy EVO 613 Li, equipped with lithium batteries, embodies the philosophy of constant improvement and innovation shared at Pride Mobility®. We believe in pushing boundaries to ensure you never fall behind.



Compact footplate



13" solid wheels



Lithium ion battery

## Specifications

Weight capacity	136 kg
Speed	Up to 8,2 km/h
Range per charge	Up to 33,8 km
Turning radius	559 mm
Ground clearance	63 mm
Length	889 - 1041 mm
Width	559 mm
Drive wheels	13" solid
Caster wheels	5" solid
Weight without batteries	64,5 kg
Batteries	Type: Lithium ion
More?	Download our specifications sheet on the website

## Features

- Mid-wheel drive technology
- Active-Trac suspension
- 13" solid drive tires
- Lithium ion battery
- Up to 8,2 km/h

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## Available colors

