



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl


The Victory Luna 3-Wheel mobility scooter delivers powerful performance both indoors and outdoors.



13" pneumatic tires



Easy to operate tiller display



Front basket standard

Specifications

Weight capacity	159 kg
Speed	Up to 15 km/h
Range per charge	Up to 41 km
Turning radius	1257 mm
Ground clearance	95 mm
Length	1397 mm
Width	699 mm
Tires	13" pneumatic, front and rear
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Weight without batteries	87,5 kg
Batteries	Type: Two 12-volt, deep-cycle Size: 55 Ah or 75 Ah
More?	Download our specifications sheet on the website

Features

- Easy throttle control by switchable push and pull lever
- Large 13" terrain gripping pneumatic tires for excellent outdoor performance
- Full lighting package including directional signals and angle-adjustable, LED head and tail lights
- Front basket standard
- Tiller display operable with just one hand

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



Red