

VICTORY *XL 130*



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl


Luxury details, such as a continuous delta handlebar and bright LED lighting, give the Victory XL 130 a streamlined and sporty design.



13" pneumatic tires



High quality design



Front and rear suspension

Specifications

Weight capacity	158,8 kg
Speed	Up to 15 km/h
Range per charge	Up to 37,9 km
Turning radius	1606 mm
Ground clearance	94 mm
Length	1355 mm
Width	635 mm
Tires	13" pneumatic, front and rear
Drive	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Weight without batteries	96 kg
Batteries	Type: Two 12-volt, sealed lead-acid, non-spillable deep-cycle Size: 55 Ah or 75 Ah
More?	Download our specifications sheet on the website

Features

- The high-visibility automatic rear brake lighting system
- Front and rear suspension
- Exclusive feather-touch disassembly provides the ultimate quick transport system
- Long-lasting, high-intensity LED headlight, taillight and indicators
- Easy access, tiller-mounted charger port for standard dual voltage off-board charger

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors

