

GO lite



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

The all-new i-Go® Lite weighs only 18 kg making it our lightest power chair yet! With a built-in USB charger, you'll be able to travel and stay connected with the world around you.



USB charger



Efficient joystick



Easy to fold

Specifications

Weight capacity	135 kg
Speed	Up to 6 km/h
Range per charge	Up to 15 km
Turning radius	900 mm
Length	940 mm
Width	610 mm
Tires	7" solid front, 8,5" solid rear
Weight without batteries	18 kg
Batteries	Type: Lithium ion Size: 12 Ah
More?	Download our specifications sheet on the website

Features

- Large secure under-seat storage compartment
- On-board USB charger
- Stylish Carbon Fibre design
- Simple one-touch folding system
- Adjustable control position

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



Carbon