




PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

Experience the beauty of lightweight travel with the i-Go®+. This compact power chair folds in just a few simple steps for effortless transport. The i-Go®+ makes traveling a breeze with its convenient storage.



Comfort headrest



Swing away controller



Underseat storage

Specifications

Weight capacity	113 kg
Speed	Up to 5,8 km/h
Range per charge	Up to 15 km
Turning radius	616 mm
Ground clearance	79 mm
Length	1016 mm
Width	597 mm
Drive wheels	12" pneumatic
Caster wheels	8" solid
Drive	Two motor, rear-wheel drive
Weight without batteries	24 kg
Batteries	Type: Lithium ion Size: 18.2 Ah
More?	Download our specifications sheet on the website

Features

- A compact and lightweight design makes traveling a breeze
- Back seat pocket for extra storage while on-the-go
- Folds in just a few simple steps for effortless transport
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extra comfort

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



Silver

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Pride Mobility Products Europe
De Zwaan 3
1601 MS Enkhuizen
The Netherlands

Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
pridemobility.eu · info@prideeurope.nl



EU_SSP0004/RevA/01MAR24