

Economy class price First class features




GO GO[®]
ULTRA X

The **Go Go Ultra X** is loaded with features such as one-hand disassembly and a convenient drop-in battery box for worry-free travel. With a 117 kg weight capacity, a maximum speed up to 6,4 km/h, a per-charge range up to 17,5 km on the 3-wheel and 18,4 km on the 4-wheel, the Go Go Ultra X is an exceptional value.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Specifications

Product number	3W: SC40X 4W: S49EU
Weight capacity	118 kg
Max. speed ¹	Up to 6,4 km/h
Ground clearance ²	32 mm
Turning radius ²	3W: 838 mm 4W: 1120 mm
Length ²	3W: 940 mm 4W: 1003 mm
Width ²	495 mm
Front tires	3W: 8" solid 4W: 7" solid
Rear tires	8" solid
Max. range per charge ^{1,4}	3W: Up to 11 km 4W: Up to 11,5 km
Total weight without batteries	3W: 34,5 kg 4W: 36,5 kg
Standard seat ²	Type: Folded Material: Black Vinyl Width: 432 mm Depth: 432 mm Weight: 10,5 kg
Drive system	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Batteries	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 12 Ah
Battery weight ⁵	8,8 kg
Color	 Red

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Features

- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Pride's exclusive black, non-scuffing tires
- Standard front basket

Accessories

- Black weather cover
- Rearview mirror
- Cup holder
- Lap belt
- Safety flag
- Saddle bag
- Cell phone holder
- XLR USB charger
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane holder, single or double