

Outdoor Adventures await



Front lights





Tiller adjustment lever

Rugged and built to last; the **Ranger** is everything you want in an outdoor scooter! Equipped with front and rear suspension, this aggressive scooter is designed for cruising on walking trails or driving around your property. Enjoy total visibility with complete LED lighting. With a user-friendly console and dual hydraulic brakes for added safety, the Ranger is the perfect way to embrace life in the great outdoors!

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



Specifications

Product number	MV600
Weight capacity	159 kg
Maximum speed ¹	Up to 15 km/h
Ground clearance ²	140 mm
Turning radius ²	2407 mm
Length ²	1556 mm
Width ²	829 mm
Tires	14,5" pneumatic
Suspension	Front and rear (CTS suspension)
Range per charge ^{1,4}	Up to 32,1 km at 159 kg Up to 38,4 km at 100 kg
Total weight without batteries ⁶	124,25 kg
Standard seating	Type: Pillow-top, limited recline Material: Black vinyl Width: 457 mm Depth: 419 mm Weight: 22 kg
Motor output	24-volt DC motor, rear-wheel drive, dual in-line
Battery requirements ^{3,6}	Quantity: 2 Type: Sealed lead-acid, deep cycle Size: 75 or 100 Ah
Battery weight	75 Ah: 24,4 kg 100 Ah: 26 kg
Battery charger	5 or 8 Amp, off-board
Colors	  Dessert Storm Stone Grey

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.
3. AGM or gel cell type required.
4. Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Incl. standard seat
6. Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

Features

- User friendly LED console displays time, temperature, miles driven and trips taken
- Adjustable delta tiller with ergonomic, wraparound handles
- USB charger built into the tiller for convenient charging of smart phones and portable devices
- Durable and stylish rear bumper
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Full LED lighting package includes headlights, hazard lights, rear backup sensor LED lights and directional signals
- XLR charging port on tiller
- Hydraulic hand brake

Accessories

- Black weather cover
- Rearview mirror
- Cup holder
- Lap belt
- Safety flag
- Saddle bag
- Cell phone holder
- XLR USB charger
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane holder, single or double

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tire pressure, ambient temperature, and the gradient and surface of the road or pavement.