

PURSUIT²



Colours Available:



Robust Features With Supreme Quality

The Colt Pursuit drives at speeds of up to 8 mph. As well as being equipped with a regenerative braking system, a comfortable high-back seat and complete LED lighting

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1460 mm (57.48")
Overall Width	753 mm (29.65")
Drive Tyres	Front: 330 mm (13") pneumatic Rear: 330 mm (13") pneumatic
Ground Clearance	96mm (3.78")
Mobility	
Turning Radius	1760 mm (57.48")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	100 mm
Dual Breaking System	Regenerative and electromechanical
Battery	
Battery Requirements	1 x 50 Ah Li 2 x 50 Ah Li (optional)
Battery Weight	10.8 kg each
Battery Charger	8 amp off-board charger
Total Weight Without Batteries	100 kg

Features



Front and rear suspension with black rims



Full lighting package



Ultra-comfort seat



REV 01 - 15/06/2023

