



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

With a mid-height backrest and pneumatic alloy wheels, the Colt Plus still offers excellent features at an affordable price.



Sleek design



10" pneumatic tires



Easy to use cockpit

Specifications

Weight capacity	159 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 25 km
Turning radius	1370 mm
Ground clearance	57 mm
Length	1194 mm
Width	560 mm
Tires	10" pneumatic, front and rear
Drive	Rear-wheel drive, single motor, 24-volt DC motor
Weight without batteries	64,5 kg
Batteries	Type: Two 12-volt, sealed lead-acid, non-spillable, deep-cycle Size: 40 Ah
More?	Download our specifications sheet on the website

Features

- Stylish 10" front and rear tires.
- The LED delta handlebar design is not only stylish, but also easy to use.
- Comfortable and practical, the mid-height backrest folds flat for easier storage and transport.

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



Red