

GO GO ULTRA



PRIDE[®]
MOBILITY

Tel: +31 (0)85-0408630 . Fax: +31 (0)85-0430039
pridemobility.eu . info@prideeurope.nl

Experience hassle-free journeys with the Go Go Ultra travel mobility scooter, designed for lightweight portability, comfort, and convenience.



8" solid tires



New delta tiller



Future proof design

Specifications

Weight capacity	120 kg
Speed	Up to 6,4 km/h
Range per charge	Up to 18,5 km
Turning radius	1385 mm
Ground clearance	32 mm
Length	1022 mm
Width	490 mm
Tires	8" solid, front and rear
Weight without batteries	36,3 kg
Batteries	Type: Lead-acid Size: 2x 12 Ah
More?	Download our specifications sheet on the website

Features

- The scooter has an ultra-lightweight and portable design, making it easy to carry and transport.
- It can be disassembled quickly, which makes storage and travel very convenient.
- The ergonomic delta tiller provides comfortable and easy control.
- There is extra legroom to ensure added comfort during longer rides.
- The new delta tiller design improves handling and ensures a more user-friendly driving experience.

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
3. Battery weight may vary based on battery manufacturer.
4. Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Available colors



PRIDE
MOBILITY

Pride Mobility Products Europe
De Zwaan 3
1601 MS Enkhuizen
The Netherlands
Tel: +31 (0)85-0408630 · Fax: +31 (0)85-0430039
pridemobility.eu · info@prideeurope.nl



EU_SSP0062/Rev0/20AUG25